# Udder Seal Dry Cow Intramammary Teat Sealant Device

# **Description:**

Udder Seal is a sterile, antibiotic-free, off-white smooth paste. It is presented, ready to use, in syringes, each containing 4 g of intramammary infusion composed of 65% Bismuth Subnitrate, in a mineral oil vehicle. Udder Seal is indicated as an aid in the prevention of new intramammary infection throughout the dry period by providing a malleable barrier in the teat canal. This results in a reduction in the incidence of clinical mastitis in the dry period.

#### For Use On:

Cattle

#### **Benefits:**

- Quick delivery
- · Accurate placement and an effective seal
- · Delivers easy syringe-ability and excellent in-teat performance
- Antibiotic-free

## Dosage and Administration:

After last milking at dry-off or during last two months before first calving in non-lactating first calf heifers, clean and disinfect the teats with the alcohol wipes. Infuse the entire contents of one syringe of Udder Seal into each quarter by inserting the nozzle into the teat and applying gentle continuous pressure to the plunger until the paste is expressed.

Do not massage teat or udder following infusion. If an antibiotic is infused into the teat, the teat should be re-wiped with an alcohol wipe prior to infusion with Udder Seal. A clean, disinfected and dry teat is essential, as Udder Seal possesses no antimicrobial activity.



### Internal teat sealants perform the function of the keratin plug

- · They form an effective seal
- They persist throughout the dry period
- · They strip out easily
- · They are insoluble to milk and harmless to calves

#### Packaging:

144 syringes per pail, 1 pail/case UPC # 7-45801-05445-5



030623

**Udder Seal Dry Cow Intramammary Teat Sealant Device** 

Tel 800-821-5570 | Fax 816-224-3080 | info@durvet.com